

Intensive whitening from 20°C. For whites that look new for longer! Dr. Beckmann Super Whitening sachets help restore that 'as new' brilliance to your whites that is just not possible using detergents alone. Dr. Beckmann's unique formula gives whites over twice as much 'glo' than using detergent alone. pH 7.5 – 8.5

#### APPLICATIONS

The last resort to achieve super white laundry.

#### WASH TEMPERATURE

20°C to 95°C



**Product** LSRDR110  
**Size** x 4 sachets

#### INSTRUCTIONS

##### Step 1

- Always check the garment's care label and machine instructions before you wash
- Wash matching items together

##### Step 2

- Place a sachet into the back of the machine drum, on top of the laundry
- Add your wash detergent to the detergent drawer as per usual
- Set programme and run machine as per usual

#### TIPS

- ⚠ Always remove empty sachet at the end of the wash cycle
- ⚠ Keep out of the reach of children

#### TO ORDER

**Shop online** [www.eyservices.com](http://www.eyservices.com)

**Call** +33 (0)4 89 02 43 40

**Email** [info@eyservices.com](mailto:info@eyservices.com)

Whilst every effort has been made to ensure the accuracy of the information supplied herein, Earthly Supply Company cannot be held responsible or liable for any damage or injury arising from the use or misuse of these products.